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CE/2010

Validity

This supplementary information describes important variations to the standard product and is only valid together with the prevailing operating instructions.

- SplitFlow 310
 - Pfeiffer Vacuum product number PM P04 037 A

This present version of a turbomolecular drag pump has been designed for a special application on the basis of the SplitFlow™ principle. This turbopump can only be operated in conjunction with its designated components.

Applicable documents

SplitFlow 310	Operating instructions
Operating instructions standard pump "HiPace 300"	PT 0202 BN*
Operating instructions "Electronic drive unit TC 110"	PT 0204 BN*
Declaration of conformity	Part of this document

*also available via www.pfeiffer-vacuum.net

Installation

DANGER

Danger from the turbopump being torn-off

In case of sudden blocking of the rotor, torques of up to 2730 Nm can occur, which can lead with incorrect attachment to tearing the turbopump off. The energy released thereby can hurl the entire pump or fragments from their inside around the area. This can cause severest injuries (possibly resulting in death) and large property damage.

- Precisely follow installation instructions.
- Only use Pfeiffer Vacuum original components (accessories) for installation.

Technical data

Parameter	SplitFlow 310
Flange (in)	
• High vacuum connection "H"	Ø 102 mm
• Intermediate stage "H1"	Ø 102 mm
• Intermediate stage "H2"	DN 25 ISO-KF
Flange (out) "VV"	DN 16 ISO-KF
Rotation speed $\pm 2\%$	60000 rpm (1000 Hz)
Standby rotation speed	40000 rpm (667 Hz)
Run-up time	10 min
Pumping speed for N ₂ on "H"	200 l/s
Pumping speed for He on "H"	210 l/s
Pumping speed for N ₂ on "H1"	250 l/s
Pumping speed for He on "H1"	220 l/s
Pumping speed for N ₂ on "H2"	20 l/s
Pumping speed for He on "H2"	22 l/s
Gas throughput at full rotation speed for N ₂ on "H"	4 mbar l/s
Gas throughput at full rotation speed for He on "H"	6 mbar l/s
Gas throughput at full rotation speed for N ₂ on "H1"	4 mbar l/s
Gas throughput at full rotation speed for He on "H1"	6 mbar l/s
Gas throughput at full rotation speed for N ₂ on "H2"	7 mbar l/s
Gas throughput at full rotation speed for He on "H2"	8 mbar l/s
Compression ratio for N ₂ (H1/H)	$2.7 \cdot 10^6$
Compression ratio for N ₂ (H2/H1)	$3.5 \cdot 10^3$
Compression ratio for N ₂ (VV/H)	$>1 \cdot 10^9$
Compression ratio for He (H1/H)	$4 \cdot 10^2$
Compression ratio for He (H2/H1)	$1.5 \cdot 10^2$
Compression ratio for He (VV/H)	$2 \cdot 10^8$
Fore-vacuum max. for N ₂ on "VV"	10 mbar

Parameter	SplitFlow 310
Weight	13.5 kg

Dimensions


Fig. 1: SplitFlow 310, PM P04 037 A

Declaration of conformity

according to the EC directive:

- **Machinery 2006/42/EC (Annex II, no. 1 A)**

We hereby declare that the product cited below satisfies all relevant provisions of EC directive "Machinery" **2006/42/EC**.

In addition, the product cited below satisfies all relevant provisions of EC directive "Electromagnetic Compatibility" **2004/108/EC**.

The agent responsible for compiling the technical documentation is Mr. Jörg Stanzel, Pfeiffer Vacuum GmbH, Berliner Straße 43, 35614 Asslar.

SplitFlow 310

Guidelines, harmonised standards and national standards and specifications which have been applied:

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Signatures:

Pfeiffer Vacuum GmbH
Berliner Straße 43
35614 Asslar
Germany

(M. Bender)
Managing Director

(Dr. M. Wiemer)
Managing Director

CE/2010

Vacuum is nothing, but everything to us!

Turbopumps

Rotary vane pumps

Roots pumps

Dry compressing pumps

Leak detectors

Valves

Components and feedthroughs

Vacuum measurement

Gas analysis

System engineering

Service

PFEIFFER  **VACUUM**

Pfeiffer Vacuum Technology AG · Headquarters/Germany

Tel. +49-(0) 64 41-8 02-0 · Fax +49-(0) 64 41-8 02-2 02 · info@pfeiffer-vacuum.de · www.pfeiffer-vacuum.net